

Intelicode® Unattended Install

To perform an Intelicode® unattended install you will need to download a copy of the latest installer for Intelicode®. This can be obtained from this page: <https://intelicode.com/update>

Next you will need to open a command prompt and execute the following:

```
Intelicode_<version number>_Setup.exe /qn
```

An example of this would be `Intelicode_17.5.0.6_Setup.exe /qn` but the version number will change with each update that is released so be sure to use the correct filename from the file you downloaded. The current directory of your command prompt will also need to be the folder your download was saved to otherwise you will need to specify the path to it as well.

The `/qn` switch sets the installer to run in quiet mode with no UI.

The following modes will also work if needed:

`/passive` - unattended mode (the installation shows only a progress bar)

`/qr` - reduced UI

A restart should not be required after installing Intelicode®.

A restart may be required for Microsoft SQL Server if you installed that during the same session. Windows should notify you is necessary.

If this process is used for an upgrade then Intelicode® would need to be closed before performing the upgrade.